

American
National
Standard

ANSI/AAMI
PC76:2021

Active implantable medical
devices—Requirements and
test protocols for safety of
patients with pacemakers
and ICDs exposed to
magnetic resonance imaging

Active implantable medical devices—Requirements and test protocols for safety of patients with pacemakers and ICDs exposed to magnetic resonance imaging

Approved 11 February 2021 by
AAMI

Approved 06 April 2021 by
American National Standards Institute, Inc.

Abstract: Provides requirements and test protocols for implantable pacemakers and ICDs exposed to magnetic resonance imaging. Physicians are increasingly using magnetic resonance imaging as tool for differential diagnostic, thus exposing pacemakers and ICD patients to such equipment. Current product standards for implantable pacemakers and ICDs do not include requirements and test protocols for implantable pacemakers and ICDs, which would ensure patient safety during such procedures.

Keywords: implantable medical devices, testing, pacemakers, ICDs, magnetic resonance, imaging

AAMI Standard

This Association for the Advancement of Medical Instrumentation (AAMI) standard implies a consensus of those substantially concerned with its scope and provisions. The existence of an AAMI standard does not in any respect preclude anyone, whether they have approved the standard or not, from manufacturing, marketing, purchasing, or using products, processes, or procedures not conforming to the standard. AAMI standards are subject to periodic review, and users are cautioned to obtain the latest editions.

CAUTION NOTICE: This AAMI standard may be revised or withdrawn at any time. AAMI procedures require that action be taken to reaffirm, revise, or withdraw this standard no later than five years from the date of publication. Interested parties may obtain current information on all AAMI standards by calling or writing AAMI.

All AAMI standards, recommended practices, technical information reports, and other types of technical documents developed by AAMI are voluntary, and their application is solely within the discretion and professional judgment of the user of the document. Occasionally, voluntary technical documents are adopted by government regulatory agencies or procurement authorities, in which case the adopting agency is responsible for enforcement of its rules and regulations.

Published by
AAMI
901 N. Glebe Road, Suite. 300
Arlington, VA 22203-1853
www.aami.org

© 2021 by the Association for the Advancement of Medical Instrumentation

All Rights Reserved

This publication is subject to copyright claims of ISO and AAMI. No part of this publication may be reproduced or distributed in any form, including an electronic retrieval system, without the prior written permission of AAMI. All requests pertaining to this draft should be submitted to AAMI. It is illegal under federal law (17 U.S.C. § 101, et seq.) to make copies of all or any part of this document (whether internally or externally) without the prior written permission of the Association for the Advancement of Medical Instrumentation. Violators risk legal action, including civil and criminal penalties, and damages of \$100,000 per offense. For permission regarding the use of all or any part of this document, contact AAMI, 901 N. Glebe Road, Suite 300, Arlington, VA 22203 Phone: (703) 525-4890; Fax: (703) 525-1421

Printed in the United States of America

ISBN 978-1-57020-813-3

Contents

Page

Committee representation	viii
Foreword	ix
1 Scope	1
2 Normative references	1
3 Terms and definitions	2
4 Symbols and abbreviated terms	6
5 General requirements for non-implantable parts	6
6 Global malfunction	6
7 General considerations for application of the tests of this document	9
8 Protection from harm to the patient from lead electrode heating	10
9 Protection from harm to the patient from device heating	25
10 Protection from harm to the patient caused by gradient-induced vibration	29
11 Protection from harm to the patient caused by B_0 -induced force	32
12 Protection from harm to the patient caused by B_0 -induced torque	37
13 Protection from harm to the patient caused by gradient-induced extrinsic electric potential	38
14 Protection from harm to the patient caused by B_0 -induced malfunction	66
15 Protection from harm to the patient caused by RF-induced malfunction and RF rectification	70
16 Protection from harm to the patient caused by gradient-induced malfunction	79
17 Combined fields test	90
Annex A (informative) RF injection method utilizing standard leads	97
Annex B (informative) Use of standard production leads to deliver a target RF power at 64 MHz to the myocardium in an animal test	100
Annex C (informative) CEM43 °C (cumulative equivalent minutes at 43 °C)	104
Annex D (informative) <i>In vivo</i> temperature rise in response to applied heat flux	108
Annex E (normative) Gradient RMS (GCS) calculation method	110
Annex F (normative) RF GCS compliance criteria	113
Annex G (informative) Example of implementation of MR scan protocols for combined fields	117
Bibliography	121
Figures	
Figure 1—Diagram showing the interval between two valid pacing amplitudes in different situations	9
Figure 2—Example E_{tan} magnitude and phase	20
Figure 3— Δ PCT versus power	22
Figure 4—Gradient heating waveform with a peak value of 127 T/s, RMS value of 42 T/s, T1= 0.2 ms, T3= 3.2 ms to be used for gradient can heating assessment	26

Figure 5—Example of thermally insulating container	27
Figure 6—Example setup with resistive heaters placed on one side of a device	28
Figure 7—General test flow for gradient-induced vibrations	30
Figure 8—MR system region of high dB/dz	33
Figure 9—Lead measurement configuration	34
Figure 10—Device surface areas	35
Figure 11—Example gradient induced PG voltage and current	39
Figure 12—Gradient test pulse synchronization to pace pulse	40
Figure 13—Gradient test pulse synchronization to pace charge balance pulse	41
Figure 14—Gradient test pulse synchronization to pace pulse and charge balance pulse	42
Figure 15—Test setup to measure gradient induced PG charge, current injection, rectification, and distortion.....	44
Figure 16—Simulated gradient-induced voltage pulse waveform	45
Figure 17—Combined gradient induced PG charge measurement.....	47
Figure 18—Separate gradient induced PG charge and current measurement	50
Figure 19—Example simulated gradient-induced voltage pulse train for rectification test.....	51
Figure 20—PG Gradient pulse train rectified current measurement.....	52
Figure 21—PG gradient pace distortion synchronization	55
Figure 22—Single chamber pacemaker test connections.....	59
Figure 23—Dual chamber pacemaker test connection.....	60
Figure 24—CRT pacemaker test connections.....	61
Figure 25—Single chamber, single shock coil ICD test connections.....	62
Figure 26—Dual chamber, two shock coil ICD test connections.....	63
Figure 27—CRT, Two shock coil ICD test connections.....	64
Figure 28—Tissue interface, conducted immunity test.....	65
Figure 29—Example of RF rectification pulse along with pacing pulse.....	71
Figure 30—Incident electric field (a.) and magnetic field (b.) in the phantom without the device present for HPM phantom.....	73
Figure 31—Incident electric field (a.) and magnetic field (b.) in the phantom without the device present for LHM phantom.....	74
Figure 32—Schematic diagram of an example test setup to evaluate RF-induced malfunction and unintended stimulation.....	75
Figure 33—Characteristics of RF rectangular modulation pulse	77
Figure 34—Characteristics of RF sinc modulation pulse.....	77

Figure 35—RF Test condition 1 pulse synchronization to pacing therapy.....	78
Figure 36—Radiated magnetic field test setup for gradient induced malfunction.....	81
Figure 37—Radiated gradient test signal.....	82
Figure 38—Injected gradient test signal.....	86
Figure 39—Phantom position on the MR scanner table within the scanner bore.....	91
Figure 40—DUT fixation inside ASTM phantom for combined field testing.....	94
Figure A.1—Temperature spatial distribution comparison between injected 500 kHz and radiated 64 MHz power.....	98
Figure A.2—Lead body and Tip/Ring in air.....	98
Figure A.3—Lead body in HEC, tip/ring in air.....	98
Figure B.1—One configuration of the RF power injection system used for validation tests and animal tests.....	102
Figure C.1—Comparison of measured chest skin surface temperatures from various sources.....	105
Figure C.2—Probability of IPG pocket starting temperature.....	106
Figure D.1—Equation estimating temperature rise (degrees Celsius) in response to applied heat flux (mW/cm ²).....	109
Figure E.1—Flowchart for assessing gradient induced charge and current limits.....	110
Figure E.2—Tissue stimulation strength duration curve.....	111
Figure E.3—Example UCS probability analysis.....	112
Figure E.4—Example determination of gradient induced charge and current limits.....	112
Figure F.1—Stress-strength analysis.....	113
Figure F.2—Maximum rectified voltage envelope.....	115
Figure F.3—Stimulation strength-duration curve.....	116

Tables

Table 1—Malfunction test and corresponding clause, and MR fields causing the assessed hazard.....	6
Table 2—Assessments to be performed for both injected and radiated malfunction tests.....	7
Table 3—RF coil geometry for in-vivo simulations.....	11
Table 4—Virtual body models.....	12
Table 5—RF coil constraints.....	17
Table 6—Acceptance criteria for pocket tissue.....	29
Table 7—Gradient pulse test signal parameters.....	46
Table 7—Maximum and minimum lead electrode voltages.....	56
Table 8—Gradient coil geometry for <i>in vivo</i> simulations.....	57
Table 9—Pacemaker lead voltage test combinations.....	58

Table 10—ICD Lead voltage test combinations	58
Table 11—Component values for tissue interface circuit	66
Table 12—Device classification	67
Table 13—Recommended pacing impedance values for RF injection	74
Table 14—Radiated gradient test signal parameters	83
Table 15—Maximum and minimum lead loop areas	84
Table 16—Pacemaker lead loop area test combinations	84
Table 17—ICD lead loop area test combinations	85
Table 18—Injected test signal parameters	87
Table 19—Maximum and minimum lead electrode voltages	88
Table 20—Pacemaker lead voltage test combinations	89
Table 21—ICD lead voltage test combinations	89
Table 22—Locations of the MR scanner iso-center (landmark position) with respect to the ASTM phantom used during the combined fields test	92
Table 23—List of scan protocols, with corresponding landmark position	92
Table 24—Fixation setup of system under test	94
Table A.1—Input vs. output power deliver with lead body in air and PEC	99
Table A.2—Current measurements in different test medium	99
Table C.1—Summary of sources for skin surface temperature in the chest region	105
Table C.2—Summary of additional references for skin surface temperature	105
Table C.3—Literature values of subcutaneous temperatures in the thorax region	106
Table C.4—Summary of CEM43°C acceptance criteria for pocket tissue	107
Table C.5—Summary of CEM43°C acceptance criteria for tissue surrounding lead body	107
Table D.1—Data used to determine functional relationship for temperature rise in response to applied heat flux	108
Table G.1—List of abbreviations used	117
Table G.2—List of MR sequences used, and corresponding MR acronym by vendor ^{a)}	118
Table G.3—Brain imaging protocol example [50] – Minimum scan time: 6 minutes	118
Table G.4—Lumbar spine imaging protocol example [50] – Minimum scan time: 4 minutes	119
Table G.5—Liver imaging protocol example [52] – Minimum scan time: 4 minutes	119
Table G.6—Kidney imaging protocol example [52] – Minimum scan time: 2 minutes	119
Table G.7—MRCP (magnetic resonance cholangiopancreatography) imaging protocol example [52] – Minimum scan time: 2 minutes	119
Table G.8—Enterography imaging protocol example [52] – Minimum scan time: 2 minutes	119

Table G.9—Breast imaging protocol example [51] – Minimum scan time: 2 minutes.....	119
Table G.10—Cardiac imaging protocol example [53] – Minimum scan time: 5 minutes.....	120
Table G.11—Hand imaging protocol example [54] – Minimum scan time: 3 minutes	120
Table G.12—Shoulder imaging protocol example [55] – Minimum scan time: 3 minutes.....	120
Table G.13—Prostate imaging protocol example [57] – Minimum scan time: 3 minutes.....	120
Table G.14—Vulva and vagina imaging protocol example [56] – Minimum scan time: 5 minutes	120

Currently in preview, click buy full version

Committee representation

Association for the Advancement of Medical Instrumentation

AAMI Pacemaker & MRI Compatibility Working Group

This AAMI American National Standard (ANS) was developed and approved by the AAMI Pacemaker & ICD MRI Compatibility Working Group.

At the time this document was published, the **AAMI Pacemaker & MRI Compatibility Working Group** had the following members:

Cochairs: Mitchell Shein
Shiloh Sison

Members: Jonathan Edmonson, Medtronic Inc Campus
Mitchell Shein, CPR Consulting, LLC
Shiloh Sison, Abbott Laboratories
Paul Stadnik, Micro Systems Engineering Inc
Bob Stevenson, Integer
Johan van den Brink, Philips
Ingmar Viohl, FDA/CDRH

Alternates: Federico DeMula, Integer
Krishna Kurpad, Micro Systems Engineering Inc
Gabriel Mouchawar, Abbott Laboratories
Sunder Rajan, FDA/CDRH
Ronald Reitan, Boston Scientific Corporation

NOTE Participation by federal agency representatives in the development of this standard does not constitute endorsement by the federal government or any of its agencies.

Foreword

The following verbal forms are used within AAMI documents to distinguish requirements from other types of provisions in the document:

- “shall” and “shall not” are used to express requirements;
- “should” and “should not” are used to express recommendations;
- “may” and “may not” are used to express permission;
- “can” and “cannot” are used as statements of possibility or capability;
- “might” and “might not” are used to express possibility;
- “must” is used for external constraints or obligations defined outside the document; “must” is not an alternative for “shall.”

Suggestions for improving this recommended practice are invited. Comments and suggested revisions should be sent to Standards, AAMI, 901 N Glebe Road, Suite 300, Arlington, VA 22203 or standards@aami.org.

NOTE This foreword does not contain provisions of the AAMI PC76, *Active implantable medical devices - Requirements and test protocols for safety of patients with pacemakers and ICDs exposed to magnetic resonance imaging* (AAMI PC76:2021), but it does provide important information about the development and intended use of the document.

Active implantable medical devices—Requirements and test protocols for safety of patients with pacemakers and ICDs exposed to magnetic resonance imaging

1 Scope

This document is applicable to transvenous pacemaker, ICD, and CRT systems intended to be used in patients who undergo a magnetic resonance scan in 1.5 T, cylindrical (circular or elliptical cross-section) bore, whole body magnetic resonance (MR) scanners operating at approximately 64 MHz with whole body coil excitation.

The tests that are specified in this document characterize interactions with the magnetic and electromagnetic fields associated with an MR scanner. The tests can be used to demonstrate device operation according to its MR Conditional labelling. The tests are not intended to be used for the routine testing of manufactured products. Some of the tests are type tests whereas others require sample size justification.

The scope of this document is limited to:

- systems that do not use sensing functions or that are programmed not to use sensing functions to affect therapy delivery during an MR scan;
- systems which have high voltage therapy disabled during an MR scan;
- devices which are implanted in the pectoral region of patients.

Requirements for non-implantable parts are outside the scope of this document.

The requirements of this document supplement or modify those of AAMI/ISO TIR10974.

NOTE 1 The document might provide useful guidance for newer therapies such as subcutaneous ICD systems, leadless pacemakers, and loop recorders, but their requirements are outside the scope of this document.

NOTE 2 Safety requirements for MR scanners can be found in IEC 60601-2-33.

2 Normative references

The following documents, in whole or in part, are normatively referenced in this document and are indispensable for its application. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

AAMI/ISO TIR10974, *Assessment of the safety of magnetic resonance imaging for patients with an active implantable medical device*

IEC 60601-2-33, *Medical electrical equipment — Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis*

ASTM F2052-15, *Standard Test Method for Measurement of Magnetically Induced Displacement Force on Medical Devices in the Magnetic Resonance Environment*