

# Technical Information Report

AAMI TIR33:2005

**Sterilization of health care  
products—Radiation—  
Substantiation of a selected  
sterilization dose—Method  $VD_{max}$**

---



Association for the Advancement  
of Medical Instrumentation

# **Sterilization of health care products—Radiation— Substantiation of a selected sterilization dose— Method $VD_{max}$**

Approved 21 December 2005 by  
Association for the Advancement of Medical Instrumentation

**Abstract:** This technical information report is intended to prepare the industry for extending application of the  $VD_{max}$  method contained in ANSI/AAMI/ISO 11137-2:2006, which is for use with selected doses of 15 kGy and 25 kGy, to the following additional doses: 17.5, 20, 22.5, 27.5, 30, 32.5, and 35 kGy.

**Keywords:** radiation,  $VD_{max}$ , substantiation

## AAMI Technical Information Report

A technical information report (TIR) is a publication of the Association for the Advancement of Medical Instrumentation (AAMI) Standards Board that addresses a particular aspect of medical technology.

Although the material presented in a TIR may need further evaluation by experts, releasing the information is valuable because the industry and the professions have an immediate need for it.

A TIR differs markedly from a standard or recommended practice, and readers should understand the differences between these documents.

Standards and recommended practices are subject to a formal process of committee approval, public review, and resolution of all comments. This process of consensus is supervised by the AAMI Standards Board and, in the case of American National Standards, by the American National Standards Institute.

A TIR is not subject to the same formal approval process as a standard. However, a TIR is approved for distribution by a technical committee and the AAMI Standards Board.

Another difference is that, although both standards and TIRs are periodically reviewed, a standard must be acted on—reaffirmed, revised, or withdrawn—and the action formally approved usually every 5 years, but at least every 10 years. For a TIR, AAMI consults with a technical committee about 5 years after the publication date (and periodically thereafter) for guidance on whether the document is still useful—that is, to check that the information is relevant or of historical value. If the information is not useful, the TIR is removed from circulation.

A TIR may be developed because it is more responsive to underlying safety or performance issues than a standard or recommended practice, or because achieving consensus is extremely difficult or unlikely. Unlike a standard, a TIR permits the inclusion of differing viewpoints on technical issues.

**CAUTION NOTICE:** This AAMI TIR may be revised or withdrawn at any time. Because it addresses a rapidly evolving field or technology, readers are cautioned to ensure that they have also considered information that may be more recent than this document.

All standards, recommended practices, technical information reports, and other types of technical documents developed by AAMI are *voluntary*, and their application is solely within the discretion and professional judgment of the user of the document. Occasionally, voluntary technical documents are adopted by government regulatory agencies or procurement authorities, in which case the adopting agency is responsible for enforcement of its rules and regulations.

Comments on this technical information report are invited and should be sent to AAMI, Attn: Standards Department, 1110 N. Glebe Road, Suite 220, Arlington, VA 22201-4795.

### *Published by*

Association for the Advancement of Medical Instrumentation  
1110 N. Glebe Road, Suite 220  
Arlington, VA 22201-4795

© 2006 by the Association for the Advancement of Medical Instrumentation  
All Rights Reserved

Publication, reproduction, photocopying, storage, or transmission, electronically or otherwise, of all or any part of this document without the prior written permission of the Association for the Advancement of Medical Instrumentation is strictly prohibited by law. It is illegal under federal law (17 U.S.C. § 101, *et seq.*) to make copies of all or any part of this document (whether internally or externally) without the prior written permission of the Association for the Advancement of Medical Instrumentation. Violators risk legal action, including civil and criminal penalties, and damages of \$100,000 per offense. For permission regarding the use of all or any part of this document, contact AAMI at 1110 N. Glebe Road, Suite 220, Arlington, VA 22201-4795. Phone: (703) 525-4890; Fax: (703) 525-1067.

Printed in the United States of America

ISBN 1-57020-270-2

## Contents

	Page
Glossary of equivalent standards .....	v
Committee representation .....	vii
Introduction .....	xi
<b>1</b> Scope .....	<b>1</b>
1.1 Inclusions .....	1
1.2 Exclusions .....	1
<b>2</b> Normative references .....	<b>1</b>
<b>3</b> Terms and definitions .....	<b>1</b>
<b>4</b> Definition and maintenance of product families for dose setting, dose substantiation, and sterilization dose auditing .....	<b>4</b>
<b>4.1</b> General .....	<b>4</b>
<b>4.2</b> Defining product families .....	<b>4</b>
<b>4.3</b> Designating of product to represent a product family for performance of a verification dose experiment or sterilization dose audit .....	<b>5</b>
<b>4.3.1</b> Product to represent a product family .....	<b>5</b>
<b>4.3.2</b> Master product .....	<b>6</b>
<b>4.3.3</b> Equivalent product .....	<b>6</b>
<b>4.3.4</b> Simulated product .....	<b>6</b>
<b>4.4</b> Maintaining product families .....	<b>6</b>
<b>4.4.1</b> Periodic review .....	<b>6</b>
<b>4.4.2</b> Modification to product and manufacturing process .....	<b>6</b>
<b>4.4.3</b> Records .....	<b>6</b>
<b>4.5</b> Effect of failure of establishment of sterilization dose or of a sterilization dose audit on a product family .....	<b>6</b>
<b>5</b> Selection and testing of product for establishing and verifying the sterilization dose .....	<b>6</b>
<b>5.1</b> Nature of product .....	<b>6</b>
<b>5.2</b> Sample item portion (SIP) .....	<b>7</b>
<b>5.3</b> Manner of sampling .....	<b>8</b>
<b>5.4</b> Microbiological testing .....	<b>8</b>
<b>5.5</b> Irradiation .....	<b>9</b>
<b>6</b> Method $VD_{max}$ —Substantiation of a selected sterilization dose .....	<b>9</b>
<b>6.1</b> Rationale .....	<b>9</b>
<b>6.2</b> Procedure for Method $VD_{max}$ for multiple production batches .....	<b>9</b>
<b>6.2.1</b> General .....	<b>9</b>
<b>6.2.2</b> Stage 1: Obtain samples of product .....	<b>9</b>
<b>6.2.3</b> Stage 2: Determine average bioburden .....	<b>10</b>
<b>6.2.4</b> Stage 3: Select sterilization dose .....	<b>10</b>
<b>6.2.5</b> Stage 4: Obtain $VD_{max}^{SD}$ .....	<b>10</b>
<b>6.2.6</b> Stage 5: Perform verification dose experiment .....	<b>11</b>
<b>6.2.7</b> Stage 6: Interpretation of results .....	<b>11</b>
<b>6.2.8</b> Confirmatory verification dose experiment .....	<b>11</b>
<b>6.3</b> Procedure for Method $VD_{max}$ for a single production batch .....	<b>12</b>
<b>6.3.1</b> Rationale .....	<b>12</b>
<b>6.3.2</b> General .....	<b>12</b>
<b>6.3.3</b> Stage 1: Obtain samples of product .....	<b>12</b>
<b>6.3.4</b> Stage 2: Determine average bioburden .....	<b>12</b>
<b>6.3.5</b> Stage 3: Selection of a sterilization dose .....	<b>13</b>
<b>6.3.6</b> Stage 4: Obtain $VD_{max}^{SD}$ .....	<b>13</b>

6.3.7	Stage 5: Perform verification dose experiment .....	14
6.3.8	Stage 6: Interpretation of results .....	14
6.3.9	Confirmatory verification dose experiment .....	14
7	Auditing sterilization dose .....	15
7.1	Purpose and frequency .....	15
7.2	Procedure for auditing a sterilization dose substantiated using Method $VD_{max}$ .....	15
7.2.1	General .....	15
7.2.2	Stage 1: Obtain samples of product .....	15
7.2.3	Stage 2: Determine average bioburden .....	15
7.2.4	Stage 3: Perform verification dose experiment .....	15
7.2.5	Stage 4: Interpretation of results .....	16
7.2.6	Confirmatory sterilization dose audit .....	16
8	Augmentation of a sterilization dose substantiated using Method $VD_{max}$ .....	17
9	Worked examples .....	17
9.1	Substantiation of a selected Dose: $SIP < 1.0$ and $SIP = 1.0$ .....	17
9.2	Dose auditing and augmentation: Worked example of a sterilization dose audit for a sterilization dose substantiated using Method $VD_{max}^{25}$ .....	19
Annex A	$VD_{max}^{SD}$ look-up tables for verification doses, SIP dose reduction factors, and dose augmentation values .....	21
	<b>Bibliography</b> .....	54
	<b>Tables</b>	
1	Nature of product items for establishing and verifying the sterilization dose .....	7
2	Examples for SIP calculation .....	8
3	Selected sterilization dose values and tabulation of verification dose values for specified ranges of average bioburden .....	10
4	Selected sterilization dose values and tabulation of verification dose values for specified ranges of average bioburden .....	13
5	Worked example of $VD_{max}^{25}$ substantiation ( $SIP < 1.0$ ) .....	18
6	$VD_{max}^{15}$ substantiation ( $SIP = 1.0$ ) .....	19
7	$VD_{max}^{25}$ dose audit (audit nonacceptance and augmentation) .....	20
A.1	$VD_{max}^{15}$ look-up table for verification doses, SIP dose reduction factors, and dose augmentation values for substantiation of a 15 kGy sterilization dose .....	21
A.2	$VD_{max}^{17.5}$ look-up table for verification doses, SIP dose reduction factors, and dose augmentation values for substantiation of a 17.5 kGy sterilization dose .....	22
A.3	$VD_{max}^{20}$ look-up table for verification doses, SIP dose reduction factors, and dose augmentation values for substantiation of a 20 kGy sterilization dose .....	23
A.4	$VD_{max}^{22.5}$ look-up table for verification doses, SIP dose reduction factors, and dose augmentation values for substantiation of a 22.5 kGy sterilization dose .....	25
A.5	$VD_{max}^{25}$ look-up table for verification doses, SIP dose reduction factors, and dose augmentation values for substantiation of a 25 kGy sterilization dose .....	28
A.6	$VD_{max}^{27}$ look-up table for verification doses, SIP dose reduction factors, and dose augmentation values for substantiation of a 27 kGy sterilization dose .....	31
A.7	$VD_{max}^{30}$ look-up table for verification doses, SIP dose reduction factors, and dose augmentation values for substantiation of a 30 kGy sterilization dose .....	35
A.8	$VD_{max}^{32.5}$ look-up table for verification doses, SIP dose reduction factors, and dose augmentation values for substantiation of a 32.5 kGy sterilization dose .....	41
A.9	$VD_{max}^{35}$ look-up table for verification doses, SIP dose reduction factors, and dose augmentation values for substantiation of a 35 kGy sterilization dose .....	47

## Glossary of Equivalent Standards

International Standards adopted in the United States may include normative references to other International Standards. For each International Standard that has been adopted by AAMI (and ANSI), the table below gives the corresponding U.S. designation and level of equivalency to the International Standard. NOTE: Documents are sorted by international designation.

Other normatively referenced International Standards may be under consideration for U.S. adoption by AAMI; therefore, this list should not be considered exhaustive.

International designation	U.S. designation	Equivalency
IEC 60601-1:2005	ANSI/AAMI ES60601-1:2005	Major technical variations
IEC 60601-1-2:2001 and Amendment 1:2004	ANSI/AAMI/IEC 60601-1-2:2001 and Amendment 1:2004	Identical
IEC 60601-2-04:2002	ANSI/AAMI DF80:2003	Major technical variations
IEC 60601-2-19:1990 and Amendment 1:1996	ANSI/AAMI I136:2004	Major technical variations
IEC 60601-2-20:1990 and Amendment 1:1996	ANSI/AAMI I151:2004	Major technical variations
IEC 60601-2-21:1994 and Amendment 1:1996	ANSI/AAMI/IEC 60601-2-21 and Amendment 1:2000 (consolidated texts)	Identical
IEC 60601-2-24:1998	ANSI/AAMI ID26:2004	Major technical variations
IEC/TR 60878:2003	ANSI/AAMI/IEC TIR60878:2003	Identical
IEC/TR 62296:2003	ANSI/AAMI/IEC TIR62296:2003	Identical
IEC 62304:2006	ANSI/AAMI/IEC 62304:2006	Identical
IEC/TR 62348:2006	ANSI/AAMI/IEC TIR62348:2006	Identical
ISO 5840:2005	ANSI/AAMI/ISO 5840:2005	Identical
ISO 7198:1998	ANSI/AAMI/ISO 7198:1998/2001/(R)2004	Identical
ISO 7199:1996	ANSI/AAMI/ISO 7199:1996/(R)2002	Identical
ISO 10993-1:2003	ANSI/AAMI/ISO 10993-1:2003	Identical
ISO 10993-2:2006	ANSI/AAMI/ISO 10993-2:2006	Identical
ISO 10993-3:2003	ANSI/AAMI/ISO 10993-3:2003	Identical
ISO 10993-4:2002 and Amendment 1:2006	ANSI/AAMI/ISO 10993-4:2002 and Amendment 1:2006	Identical
ISO 10993-5:1999	ANSI/AAMI/ISO 10993-5:1999	Identical
ISO 10993-6:1994	ANSI/AAMI/ISO 10993-6:1995/(R)2001	Identical
ISO 10993-7:1995	ANSI/AAMI/ISO 10993-7:1995/(R)2001	Identical
ISO 10993-9:1999	ANSI/AAMI/ISO 10993-9:1999/(R)2005	Identical
ISO 10993-10:2002 and Amendment 1:2006	ANSI/AAMI BE78:2002 ANSI/AAMI BE78:2002/A1:2006	Minor technical variations Identical
ISO 10993-11:2006	ANSI/AAMI/ISO 10993-11:2006	Identical
ISO 10993-12:2002	ANSI/AAMI/ISO 10993-12:2002	Identical
ISO 10993-13:1998	ANSI/AAMI/ISO 10993-13:1999/(R)2004	Identical
ISO 10993-14:2001	ANSI/AAMI/ISO 10993-14:2001	Identical
ISO 10993-15:2000	ANSI/AAMI/ISO 10993-15:2000	Identical
ISO 10993-16:1997	ANSI/AAMI/ISO 10993-16:1997/(R)2003	Identical
ISO 10993-17:2002	ANSI/AAMI/ISO 10993-17:2002	Identical
ISO 10993-18:2005	ANSI/AAMI BE83:2006	Major technical variations
ISO/TS 10993-19:2006	ANSI/AAMI/ISO TIR10993-19:2006	Identical

International designation	U.S. designation	Equivalency
ISO/TS 10993-20:2006	ANSI/AAMI/ISO TIR10993-20:2006	Identical
ISO 11135:1994	ANSI/AAMI/ISO 11135:1994	Identical
ISO 11137-1:2006	ANSI/AAMI/ISO 11137-1:2006	Identical
ISO 11137-2:2006 (2006-08-01 corrected version)	ANSI/AAMI/ISO 11137-2:2006	Identical
ISO 11137-3:2006	ANSI/AAMI/ISO 11137-3:2006	Identical
ISO 11138-1: 2006	ANSI/AAMI/ISO 11138-1:2006	Identical
ISO 11138-2: 2006	ANSI/AAMI/ISO 11138-2:2006	Identical
ISO 11138-3: 2006	ANSI/AAMI/ISO 11138-3:2006	Identical
ISO 11138-4: 2006	ANSI/AAMI/ISO 11138-4:2006	Identical
ISO 11138-5: 2006	ANSI/AAMI/ISO 11138-5:2006	Identical
ISO/TS 11139:2006	ANSI/AAMI/ISO 11139:2006	Identical
ISO 11140-1:2005	ANSI/AAMI/ISO 11140-1:2005	Identical
ISO 11140-5:2000	ANSI/AAMI ST66:1999	Major technical variations
ISO 11607-1:2006	ANSI/AAMI/ISO 11607-1:2006	Identical
ISO 11607-2:2006	ANSI/AAMI/ISO 11607-2:2006	Identical
ISO 11737-1: 2006	ANSI/AAMI/ISO 11737-1:2006	Identical
ISO 11737-2:1998	ANSI/AAMI/ISO 11737-2:1998	Identical
ISO 11737-3:2004	ANSI/AAMI/ISO 11737-3:2004	Identical
ISO 13485:2003	ANSI/AAMI/ISO 13485:2003	Identical
ISO 13488:1996	ANSI/AAMI/ISO 13488:1996	Identical
ISO 14155-1:2003	ANSI/AAMI/ISO 14155-1:2003	Identical
ISO 14155-2:2003	ANSI/AAMI/ISO 14155-2:2003	Identical
ISO 14160:1998	ANSI/AAMI/ISO 14160:1998	Identical
ISO 14161:2000	ANSI/AAMI/ISO 14161:2000	Identical
ISO 14937:2000	ANSI/AAMI/ISO 14937:2000	Identical
ISO/TR 14969:2004	ANSI/AAMI/ISO TIR14969:2004	Identical
ISO 14971:2000 and A1:2003	ANSI/AAMI/ISO 14971:2000 and A1:2003	Identical
ISO 15223:2000, A1:2002, and A2:2004	ANSI/AAMI/ISO 15223:2000, A1:2001, and A2:2004	Identical
ISO 15225:2000 and A1:2004	ANSI/AAMI/ISO 15225:2000/(R)2006 and A1:2004/(R)2006	Identical
ISO 15674:2001	ANSI/AAMI/ISO 15674:2001	Identical
ISO 15675:2001	ANSI/AAMI/ISO 15675:2001	Identical
ISO/TS 15843:2000	ANSI/AAMI/ISO TIR15843:2000	Identical
ISO 15882:2003	ANSI/AAMI/ISO 15882:2003	Identical
ISO/TR 16142:2006	ANSI/AAMI/ISO TIR16142:2006	Identical
ISO 17664:2004	ANSI/AAMI ST81:2004	Major technical variations
ISO 17665-1:2006	ANSI/AAMI/ISO 17665-1:2006	Identical
ISO 18472:2006	ANSI/AAMI/ISO 18472:2006	Identical
ISO/TS 19218:2005	ANSI/AAMI/ISO 19218:2005	Identical
ISO 25539-1:2003 and A1:2005	ANSI/AAMI/ISO 25539-1:2003 and A1:2005	Identical

<sup>1</sup>In production

<sup>2</sup>Final approval pending

## Committee representation

### Association for the Advancement of Medical Instrumentation Sterilization Standards Committee

The development of AAMI TIR33:2005 as an AAMI technical information report was initiated by the AAMI Radiation Sterilization Working Group (ST/WG 02) under the auspices of the AAMI Sterilization Standards Committee (AAMI/ST).

Committee approval of this document does not necessarily imply that all committee members voted for its approval.

At the time this document was published, the **AAMI Radiation Sterilization Working Group** had the following members:

<i>Cochairs</i>	Lisa Foster Byron J. Lambert
<i>Members</i>	Chris B. Anderson, St. Jude Medical Inc. Leonard S. Berman, PhD, Pall Medical Paul William Boentges, BS, Cardinal Health (MP&S) Anne F. Booth, MS, Conmed Corp. John Broad, NAMSA Trabue D. Bryans, AppTec Virginia C. Chamberlain, (Independent Expert) Rod Chu, MDS Nordion Gary N. Cranston, Consulting & Technical Svcs/PCS Greg Crego, Ethox Corp. Kate Davenport, Northview Biosciences Douglas D. Davie, Sterilization Validation Services Brian R. Drumheller, CR Bard Barry P. Fairand, PhD, (Independent Expert) Lisa Foster, Sterigenics International Ruth Garcia, Steris Corporation Joyce M. Hansen, JM Hansen & Associates Thomas L. Hansen, Terumo Medical Corporation Doug F. Harbrecht, Boston Scientific Corporation Deborah A. Havlik, Hospira Inc. Arthur H. Heiss, PhD, Bueker Plastics Corp. Craig M. Herring, Johnson & Johnson Carolyn L. Kinsley, Pharmaceutical Systems Inc. Byron J. Lambert, F. I. D. Liquidant Corporation/Cardiac Rhythm Management Jeff Martin, Alcor Laboratories Inc. David Ford McGarrick, BS, Abbott Laboratories James E. McGowan Jr., BSMBA, Sterile Works Inc. Joseph M. Meier, Ethide Laboratories Inc. Russell D. Mills, Zimmer Inc. Gerry A. O'Dell, MS, Gerry O'Dell Consulting Frank Peacock Jr., Bausch & Lomb Inc. Clyde Pearce, Clearant Inc. Manuel Saavedra Jr., Kimberly-Clark Corporation Genius V. Seliokas, Stericon Inc. Jon Seulean, Cobe Sterilization Services Inc. Harry L. Shaffer, Sterilization Consulting Services William N. Thompson, TYCO Healthcare/Kendall Richard L. Weisman, Fresenius Medical Care NA Dialysis Products Division Patrick B. Weixel, FDA/CDRH James L. Whitby, MA, MB, FRCP, University of Western Ontario (Independent Expert) Thelma Wilcott, Becton Dickinson & Company John Andrew Williams, BS, Baxter Healthcare Corporation Martell Kress Winters, B.S., SM, Nelson Laboratories Inc
<i>Alternates</i>	Lisa Baryschpolec, Johnson & Johnson Harry F. Bushar, PhD, FDA/CDRH

Charles Cogdill, Boston Scientific Corporation  
 John DiCaro, TYCO Healthcare/Kendall  
 Joyce Kay Elkins, Zimmer Inc  
 Niki Fidopiastis, Sterigenics International  
 Donna Horner, Guidant Corporation/Cardiac Rhythm Management  
 Joseph A. Hutson, Cardinal Health (MP&S)  
 Orlando C. Johnson, Hospira Inc  
 Bert Kingsbury, Terumo Medical Corporation  
 Ezra A. Koski, Cobe Sterilization Services Inc.  
 Mary Malarkey, FDA/CDRH  
 Mary S. Mayo, CR Bard  
 Consuelo Lorraine McChesney, BS, Alcon Laboratories Inc.  
 Patrick J. McCormick, PhD, Bausch & Lomb Inc.  
 Dave Parente, NAMSA  
 Timothy Ramsey, BS, Northview Biosciences  
 Mark Seybold, Baxter Healthcare Corporation  
 Ralph Stick, AppTec  
 John W. Walker, Steris Corporation  
 Wendy Wangsgard, PhD, Nelson Laboratories Inc.  
 David Weppner, Ethox Corp.

**AAMI Sterilization Standards Committee**

*Cochairs*

Victoria M. Hitchins  
 William E. Young

*Members*

Trabue D. Bryans, AppTec  
 Virginia C. Chamberlain, PhD, (Independent Expert)  
 Nancy Chobin, RN, CSPDM, (Independent Expert)  
 Anne M. Cofield, CRCST, FCS, IAHCSSM  
 Charles Cogdill, Boston Scientific Corporation  
 Ramona Conner, RN, MSN, CNOR, Association of Homeoperative Registered Nurses  
 Jacqueline Daley, Association for Professionals in Infection Control and Epidemiology  
 Kimbrell Darnell, CR Bard  
 Lisa Foster, Sterigenics International  
 James M. Gibson Jr., JM Gibson Associates  
 Barbara J. Goodman, RN, BS, CNOR (Independent Expert)  
 Joel R. Gorski, PhD, NAMSA  
 Deborah A. Havlik, Hospira Inc.  
 Victoria M. Hitchins, PhD, FDA/CDRH  
 Richard M. Johnson, MSc, BSc, Abbott Laboratories  
 Lois Atkinson Jones, MS (Independent Expert)  
 John Kowalski, PhD, MicroGammMA LLC  
 Byron J. Lambert, PhD, Guidant Corporation/Cardiac Rhythm Management  
 Colleen Patricia Landers, RN, Canadian Standards Association  
 David Liu, Johnson & Johnson  
 Jeff Martin, Alcon Laboratories Inc.  
 Patrick James McCormick, PhD, Bausch & Lomb Inc.  
 Thomas C. Moore, Getinge USA  
 Barry F. J. Page, (Independent Expert)  
 Nancy J. Rakiewicz, Ethox Corp.  
 Neil M. Schneider, 3M Healthcare  
 Michael H. Scholla, Dupont Nonwovens  
 Mark Seybold, Baxter Healthcare Corporation  
 Andrew Sharavara, Propper Manufacturing Co. Inc.  
 Frank Sizemore, American Society for Healthcare Central Service Professionals  
 Gregory O. Stecklein, MS MSM, Cardinal Health (MP&S)  
 William N. Thompson, TYCO Healthcare/Kendall  
 John W. Walker, Steris Corporation  
 James L. Whitby, MA, MB, FRCP, University of Western Ontario (Independent Expert)  
 Thelma Wilcott, Becton Dickinson & Company  
 Martell Kress Winters, BSSM, Nelson Laboratories Inc.  
 William E. Young (Independent Expert)

*Alternates*

Lloyd Brown, TYCO Healthcare/Kendall  
 Lina C. Bueno, Dupont Nonwovens

Craig M. Herring, Johnson & Johnson  
Clark W. Houghtling, Steris Corporation  
Danny Hutson, Cardinal Health (MP&S)  
Jim Kaiser, Bausch & Lomb Inc.  
Susan G. Klacik, AS, BS, IAHCSSM  
Joseph J. Lasich, BS, Alcon Laboratories Inc.  
Chiu Lin, PhD, FDA/CDRH  
Lisa N. Macdonald, Becton Dickinson & Company  
Ralph Makinen, Guidant Corporation/Cardiac Rhythm Management  
Mary S. Mayo, CR Bard  
David Ford McGoldrick, BS, Abbott Laboratories  
Jerry R. Nelson, MS PhD, Nelson Laboratories Inc.  
Jeff Peltier, Boston Scientific Corporation  
Janet Prust, 3M Healthcare  
Mike Sadowski, Baxter Healthcare Corporation  
Ralph Stick, AppTec  
Jason Voisinet, Ethox Corp.  
Valerie Welter, Hospira Inc.  
William T. Young, Sterigenics International

---

NOTE—Participation by federal agency representatives in the development of this document does not constitute endorsement by the federal government or any of its agencies.

---

## Introduction

This document is intended to be used in conjunction with ANSI/AAMI/ISO 11137, *Sterilization of health care products Radiation sterilization* (series). One of the activities encompassed within the standard is the selection of the sterilization dose to be applied to health care products.

ANSI/AAMI/ISO 11137-2:2006 includes Method  $VD_{max}$  for the substantiation of 25 kGy as a sterilization dose (termed  $VD_{max}^{25}$ ) for products with an average bioburden less than or equal to 1,000, and  $VD_{max}^{15}$  for the substantiation of 15 kGy as a sterilization dose for products with low average bioburden (less than or equal to 1.5).

This TIR parallels and extends ANSI/AAMI/ISO 11137-2; it provides a methodology for the substantiation of a range of sterilization doses. Included are 15, 17.5, 20, 22.5, 25, 27.5, 30, 32.5, and 35 kGy; each of these sterilization doses is valid only for a specified unique range of average bioburden on product.

The method described in this report is a simple approach to dose substantiation to attain a sterility assurance level (SAL) of  $10^{-6}$  at the selected sterilization dose, (e.g. 25 kGy for  $VD_{max}^{25}$ ). The application of this method is not limited by batch size or production frequency, and the number of product units irradiated in the verification dose experiment remains constant. The method employs as its basis the standard distribution of resistances (SDR) on which Method 1 is also founded and embodies the following three principles:

- Existence of a direct link between the outcome of the verification dose experiment and the attainment of an SAL of  $10^{-6}$  at the sterilization dose.
- Possession of a level of conservativeness at least equal to that of the SDR.
- For a given bioburden, use of a maximum verification dose ( $VD_{max}$ ) commensurate with substantiation of a selected sterilization dose.

This approach to sterilization dose substantiation was first outlined by Kowalski and Tallentire (1999), and from subsequent evaluations involving computational techniques (Kowalski, Aoshuang, and Tallentire, 2000), it was concluded that the method is safe and yields unambiguous results. An overview of the method is provided in Kowalski and Tallentire (2000); application of the  $VD_{max}$  approach to doses other than 25 kGy is discussed in Kowalski and Tallentire (2003).

The method described here and designated  $VD_{max}$  procedurally comprises elements that closely parallel those of dose-setting Method 1. One key area of difference is the number of product units used in the verification dose experiment. In the computer evaluations referred to above, changing the verification SAL value had little effect on the substantiation outcome, and this finding led to a sample size of 10 product units being chosen for subsequent field evaluations and, ultimately, as the basis of the sampling plan described in this document. Manufacturers of health care products who intend to use the protocols contained in this technical information report are reminded that the requirements for all users of radiation sterilization contained in ANSI/AAMI/ISO 11137 apply to the manufacture and control of production batches for which a sterilization dose is substantiated by this method. In particular, one requirement states that products shall be manufactured in circumstances such that the bioburden is controlled. Compliance with the requirements for properly controlling the quality of raw materials and the manufacturing environment, and for establishing the basic properties of the packaging material, is essential.

# Sterilization of health care products— Radiation—Substantiation of a selected sterilization dose—Method $VD_{max}$

## 1 Scope

### 1.1 Inclusions

This technical information report (TIR) describes a method of substantiating a selected dose for a sterility assurance level (SAL) of  $10^{-6}$  for radiation sterilization of health care products. This TIR also specifies a method of dose auditing to demonstrate the continued effectiveness of the sterilization dose.

NOTE 1—This method of sterilization dose substantiation may be used to meet the product qualification requirements specified in ANSI/AAMI/ISO 11137.

NOTE 2—This TIR is considered informative, and use of the terms shall, should, and so forth should be considered only within the context of this TIR. That is, if the decision is made to use this method of substantiation, then this method should be followed according to the requirements (shall) and recommendations (should) set forth in this TIR.

### 1.2 Exclusions

This method, as described, is for the substantiation of the specified doses only and cannot be used to substantiate other sterilization doses. The method cannot be used for the substantiation of a selected dose when the average bioburden estimate of the entire product unit does not meet the requirements for the selected dose.

## 2 Normative references

The following normative documents contain provisions that, through reference in this text, constitute provisions of this TIR. For dated references, subsequent amendments to or revisions of any of these publications do not apply. However, parties to agreements based on this TIR are encouraged to investigate the possibility of applying the most recent editions of the normative documents indicated below. For undated references, the latest edition of the normative document referred to applies. The Association for the Advancement of Medical Instrumentation maintains a register of currently valid AAMI standards.

ANSI/AAMI/ISO 11137-1, *Sterilization of health care products—Radiation—Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices*

ANSI/AAMI/ISO 11137-2, *Sterilization of health care products—Radiation—Part 2: Establishing the sterilization dose*

ANSI/AAMI/ISO 11137-3, *Sterilization of health care products—Radiation—Part 3: Guidance on dosimetric aspects*

ANSI/AAMI/ISO 11737-1, *Sterilization of medical devices—Microbiological methods—Part 1: Determination of population of microorganisms on products*

ANSI/AAMI/ISO 11737-2, *Sterilization of medical devices—Microbiological methods—Part 2: Tests of sterility performed in the validation of a sterilization process*

## 3 Terms and definitions

For the purposes of this TIR the following terms and definitions apply.

**3.1 absorbed dose:** Quantity of ionizing radiation energy imparted per unit mass of matter.

NOTE—The unit of absorbed dose is the gray (Gy), where 1 gray is equivalent to absorption of 1 joule per kilogram.

[ANSI/AAMI/ISO 11137-1]