



BSI Standards Publication

**Food analysis — Determination
of pesticide residues by
GC-MS — Retention times,
mass spectrometric parameters
and detector response
information**

NO COPYING WITHOUT BSI PERMISSION EXCEPT AS PERMITTED BY COPYRIGHT LAW

National foreword

This Published Document is the UK implementation of CEN/TR 16468:2013.

The UK participation in its preparation was entrusted to Technical Committee AW/275, Food analysis - Horizontal methods.

A list of organizations represented on this committee can be obtained on request to its secretary.

This publication does not purport to include all the necessary provisions of a contract. Users are responsible for its correct application.

© The British Standards Institution 2013.

Published by BSI Standards Limited 2013

ISBN 978 0 580 79327 1

ICS 65.100.01; 67.050

Compliance with a British Standard cannot confer immunity from legal obligations.

This Published Document was published under the authority of the Standards Policy and Strategy Committee on 30 April 2013.

Amendments issued since publication

Date	Text affected
------	---------------

ICS 65.100.01; 67.050

English Version

**Food analysis - Determination of pesticide residues by GC-MS -
Retention times, mass spectrometric parameters and detection
response information**

Analyse des produits alimentaires - Détermination des
résidus de pesticides par GC-SM - Temps de rétention,
paramètres de spectrométrie de masse et information sur
la réponse des détecteurs

Lebensmitteluntersuchung - Bestimmung von
Pestizidrückständen mit GC-MS - Retentionszeiten,
Parameter für die Massenspektrometrie und
Detektorempfindlichkeit

This Technical Report was approved by CEN on 3 December 2012. It has been drawn up by the Technical Committee CEN/TC 275.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, Former Yugoslav Republic of Macedonia, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

Management Centre: Avenue Marnix 17, B-1000 Brussels

Contents

Page

Foreword.....	3
Introduction	4
1 Scope.....	5
2 Analyte specific parameters for gas chromatographic determination of pesticides	5
2.1 General	5
2.2 GC-MS, GC-ECD and GC-NPD Parameters	5
2.3 Analyte specific retention times and MS parameters.....	6
Bibliography.....	24

Foreword

This document (CEN/TR 16468:2013) has been prepared by Technical Committee CEN/TC 275 "Food analysis - Horizontal methods", the secretariat of which is held by DIN.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN [and/or CENELEC] shall not be held responsible for identifying any or all such patent rights.

Currently in preview, click buy full version

Introduction

Pesticide residue analysis employs multiple methods involving extraction of residues from foods and clean-up of the extract to obtain as many analytes as possible in the purified extracts. Afterwards the extracts can be analysed by different kinds of instruments.

The hyphenation of gas chromatography (GC) and mass spectrometry (MS) is one of the most often used universal and selective analysis techniques for identification and quantification of pesticide residues in food extracts.

For the ionisation of the analytes (pesticides and/or their metabolites) in GC-MS, electron impact ionisation (EI) is most commonly used. If the typical electron energy of 70 eV is used, very often molecular ions (cation radicals) and several fragment ions are formed simultaneously.

The selective determination of each target analyte is performed by simultaneous acquisition of typically three ions formed by the analyte in the selected ion monitoring (SIM) mode. The diagnostic value of an ion depends on its mass. Usually even-numbered, high mass ions are less frequently formed and their recording results in more specific chromatograms. A reduction of selectivity may be caused by background ions formed during 'column bleeding' or by sample matrix (e.g. typical ions from fatty acids or hydrocarbons).

1 Scope

This Technical Report lists mass spectrometric parameters which are useful for the application of European Standards for the determination of pesticide residues in foods that use GC-MS, such as the following standards:

- EN 1528 (all parts), *Fatty food — Determination of pesticides and polychlorinated biphenyls (PCBs)*;
- EN 12393 (all parts), *Foods of plant origin — Multiresidue methods for the gas chromatographic determination of pesticide residues*;
- EN 15662, *Foods of plant origin — Determination of pesticide residues using GC-MS and/or GC-MS/MS following acetonitrile extraction/partitioning and clean-up by dispersive SPE — QuEChERS-method*.

To facilitate the determination of pesticides and/or metabolites using GC-MS, Table 1 specifies the diagnostic ions suitable for quantification, which can be used.

2 Analyte specific parameters for gas chromatographic determination of pesticides

2.1 General

All values indicated in Table 1 were acquired using GC-MS, GC-ECD or GC-NPD systems under the experimental conditions as outlined in 2.2. Comparative investigation showed that these parameters can be transferred simply on instruments of other types of the same or other manufacturers [1].

2.2 GC-MS, GC-ECD and GC-NPD Parameters

The following GC operating conditions have been proven to be satisfactory. This is an example for appropriate experimental conditions. Equivalent conditions may be used if they can be shown to lead to the same results.

Gas chromatograph	Agilent model 6890
Carrier gas	Helium, constant flow 1,0 ml/min (typical pressure at 70 °C is 61 kPa)
Injection technique	Pulsed splitless, 200 kPa; Pulse time 1,0 min
Injector temperature	40 °C
Injection liner	Single taper, splitless (Agilent 5181-3316)
Injection volume	1 µl or 2 µl
Purge Gas	Helium; purge flow to split vent 50 ml/min; Purge time 1,5 min
Column	Fused silica capillary column HP-5MS, length 30 m, inner diameter 0,25 mm, film thickness 0,25 µm (Agilent Nr. 19091S-433)
Temperature programme	2 min 70 °C, programmed to rise at 25 °C/min to 170 °C, at 3 °C/min to 210 °C, at 30 °C/min to 290 °C, then isothermal for 9 min
Transfer Line Temperature	280 °C
Mass spectrometer	Agilent MSD 5973N inert
MS temperatures	Quadrupole 150 °C; Source 230 °C;
Ionisation	Electron impact 70 eV
Solvent Delay	4,0 min
Retention time of parathion	15,85 min (parathion is used for calculation of relative retention times)