



# EP14

## Evaluation of Commutability of Processed Samples

This document provides guidance for evaluating the commutability of processed samples by determining if they behave differently than unprocessed patient samples when two quantitative measurement procedures are compared.

.....  
A guideline for global application developed through the Clinical and Laboratory Standards Institute consensus process.

# Clinical and Laboratory Standards Institute

*Setting the standard for quality in medical laboratory testing around the world.*

The Clinical and Laboratory Standards Institute (CLSI) is a not-for-profit membership organization that brings together the varied perspectives and expertise of the worldwide laboratory community for the advancement of a common cause: to foster excellence in laboratory medicine by developing and implementing medical laboratory standards and guidelines that help laboratories fulfill their responsibilities with efficiency, effectiveness, and global applicability.

## Consensus Process

Consensus—the substantial agreement by materially affected, competent, and interested parties—is core to the development of all CLSI documents. It does not always connote unanimous agreement but does mean that the participants in the development of a consensus document have considered and resolved all relevant objections and accept the resulting agreement.

## Commenting on Documents

CLSI documents undergo periodic evaluation and modification to keep pace with advances in technologies, procedures, methods, and protocols affecting the laboratory or health care.

CLSI's consensus process depends on experts who volunteer to serve as contributing authors and/or as participants in the reviewing and commenting process. At the end of each comment period, the committee that developed the document is obligated to review all comments, respond in writing to all substantive comments, and revise the draft document as appropriate.

Comments on published CLSI documents are equally essential and may be submitted by anyone, at any time, on any document. All comments are managed according to the consensus process by a committee of experts.

## Appeal Process

When it is believed that an objection has not been adequately considered and responded to, the process for appeal, documented in the *CLSI Standards Development Policies and Processes*, is followed.

All comments and responses submitted on draft and published documents are retained on file at CLSI and are available upon request.

## Get Involved—Volunteer!

Do you use CLSI documents in your workplace? Do you see room for improvement? Would you like to get involved in the revision process? Or maybe you see a need to develop a new document for an emerging technology? CLSI wants to hear from you. We are always looking for volunteers. By donating your time and talents to improve the standards that affect your own work, you will play an active role in improving public health across the globe.

For additional information on committee participation or to submit comments, contact CLSI.

Clinical and Laboratory Standards Institute

P: +1.610.688.0100

F: +1.610.688.0700

[www.clsi.org](http://www.clsi.org)

[standard@clsi.org](mailto:standard@clsi.org)

---

## Evaluation of Commutability of Processed Samples

Karl De Vore, BA, SSBB  
Yashpal Agrawal, MD, PhD  
Todd D. Alspach, MT(ASCP), MOS  
Jeffrey R. Budd, PhD  
Ramon A. Durazo-Arvizu, PhD  
John H. Eckfeldt, MD, PhD  
Kathie Goodwin, RAC, MBA, MT(ASCP)BB

Abdel-Baset Halim, PharmD, PhD, DABCC  
Thomas A. Long  
W. Gregory Miller, PhD  
Nisar Pampori, PhD  
Justin Thomas  
Jeffrey E. Vaks, PhD  
Hubert W. Vesper, PhD

### Abstract

---

Clinical and Laboratory Standards Institute document EP14—*Evaluation of Commutability of Processed Samples* was developed for manufacturers and providers of proficiency testing or external quality assessment programs, although it is useful to clinical laboratories as well. The document helps users 1) determine whether noncommutability is the source of unexpected results that are sometimes observed with processed samples when two quantitative measurement procedures are compared, 2) display the magnitude of the effects, and 3) ensure that laboratory performance is evaluated fairly if noncommutability is present. The suggested protocol was developed using unprocessed patient samples as the standard of comparison.

Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Commutability of Processed Samples*. 4th ed. CLSI guideline EP14 (ISBN 978-1-68440-156-7 [Print]; ISBN 978-1-68440-157-4 [Electronic]). Clinical and Laboratory Standards Institute, USA, 2022.

The Clinical and Laboratory Standards Institute consensus process, which is the mechanism for moving a document through two or more levels of review by the health care community, is an ongoing process. Users should expect revised editions of any given document. Because rapid changes in technology may affect the procedures, methods, and protocols in a standard or guideline, users should replace outdated editions with the current editions of CLSI documents. Current editions are listed in the CLSI catalog and posted on our website at [www.clsi.org](http://www.clsi.org).

**If you or your organization is not a member and would like to become one, or to request a copy of the catalog, contact us at:**

**P:** +1.610.688.0100 **F:** +1.610.688.0700 **E:** [customerservice@clsi.org](mailto:customerservice@clsi.org) **W:** [www.clsi.org](http://www.clsi.org)

Copyright ©2022 Clinical and Laboratory Standards Institute. Except as stated below, any reproduction of content from a CLSI copyrighted standard, guideline, companion product, or other material requires express written consent from CLSI. All rights reserved. Interested parties may send permission requests to [permissions@clsi.org](mailto:permissions@clsi.org).

CLSI hereby grants permission to each individual member or purchaser to make a single reproduction of this publication for use in its laboratory procedure manual at a single site. To request permission to use this publication in any other manner, e-mail [permissions@clsi.org](mailto:permissions@clsi.org).

## Suggested Citation

CLSI. *Evaluation of Commutability of Processed Samples*. 4th ed. CLSI guideline EP14. Clinical and Laboratory Standards Institute; 2022.

Previous Editions:

**Proposed Guideline**

April 1999

**Approved Guideline**

March 2001

**Approved Guideline—Second Edition**

January 2005

**Approved Guideline—Third Edition**

August 2014

EP14-Ed4

ISBN 978-1-3440-156-7 (Print)

ISBN 978-1-68440-157-4 (Electronic)

ISSN 1558-6502 (Print)

ISSN 2162-2914 (Electronic)

Volume 42, Number 17

## Committee Membership

### Consensus Council

The Consensus Council sets priorities for CLSI standards development and votes on Final Draft documents to confirm that process requirements have been met. Consensus Council members are listed on the CLSI website: <https://clsi.org/standards-development/consensus-council/>

### Document Development Committee on Evaluation of Matrix Effects

**Karl De Vore, BA, SSBB**  
**Chairholder**  
**Bio-Rad Laboratories, Inc.**  
**USA**

Ramon A. Durazo-Arvizu, PhD  
 Loyola University Medical Center  
 USA

Nisar Pampori, PhD  
 FDA Center for Devices and  
 Radiological Health  
 USA

**Dorothy J. Ball, PhD**  
**Vice-Chairholder**  
**USA**

Thomas A. Long  
 College of American Pathologists  
 USA

Hubert W. Vesper, PhD  
 Centers for Disease Control and  
 Prevention  
 USA

Jeffrey R. Budd, PhD  
 Beckman Coulter  
 USA

### Expert Panel on Evaluation Protocols

Expert Panel volunteers support the development of CLSI documents by providing technical expertise in specialty areas. Expert Panel members are listed by area of expertise on the CLSI website: <https://clsi.org/standards-development/expert-panels/>

### Staff

Clinical and Laboratory Standards  
 Institute  
 USA

Laura Martin  
*Editorial Manager*

Kristy L. Leirer, MS  
*Editor*

Tabitha Kern, MS, MLS(ASCP)<sup>CM</sup>  
*Project Manager*

Catherine E.M. Jenkins, ELS  
*Editor*

Lisa M.W. Walker, MS, ELS  
*Editor*

## Acknowledgment

---

CLSI, the Consensus Council, the Expert Panel on Evaluation Protocols, the Consensus Committee on Evaluation Protocols, and the Document Development Committee on Evaluation of Matrix Effects gratefully acknowledge the following volunteers for their important contributions to the development of this guideline in 2014 and its limited revision in 2022:

Yashpal Agrawal, MD, PhD  
New York Presbyterian Hospital –  
Cornell Campus  
USA

John H. Eckfeldt, MD, PhD  
University of Minnesota  
Medical Center – Fairview  
USA

W. Gregory Miller, PhD  
Virginia Commonwealth University  
Health System  
USA

Todd D. Alspach, MT(ASCP), MOS  
Centers for Disease Control and  
Prevention  
USA

Kathie Goodwin, RAC, MBA, MT(ASCP)  
BB  
Roche Diagnostics, Inc.  
USA

Justin Thomas  
Siemens Healthcare Diagnostics Inc.  
USA

Jeffrey R. Budd, PhD  
Beckman Coulter  
USA

Abdel-Baset Halim, PharmD, PhD,  
DABCC  
DrugScan  
USA

Jeffrey E. Vaks, PhD  
Roche Molecular Systems  
USA

Karl De Vore, BA, SSBB  
Bio-Rad Laboratories, Inc.  
USA

# Contents

Abstract .....	i
Committee Membership .....	iii
Foreword .....	vii
<b>Chapter 1: Introduction .....</b>	<b>1</b>
1.1 Scope .....	1
1.2 Background .....	3
1.3 Standard Precautions .....	5
1.4 Terminology .....	5
<b>Chapter 2: Commutability Assessment .....</b>	<b>9</b>
2.1 Assessment Methods .....	10
2.2 Process Flow Chart .....	12
2.3 Materials and Samples Assembly .....	13
2.4 Measurement Procedures .....	14
2.5 Data Analysis .....	14
2.6 Documenting Results of the Commutability Evaluation .....	20
<b>Chapter 3: Conclusion .....</b>	<b>21</b>
<b>Chapter 4: Supplemental Information .....</b>	<b>23</b>
<b>References</b> .....	24
<b>Additional Resources</b> .....	26
<b>Appendix A.</b> Description of Mathematical Model Used for Evaluating Commutability of Processed Samples Using Deming Regression .....	27
<b>Appendix B.</b> Outlier Evaluation for a Measurement Procedure Comparison Using Deming Regression .....	33
<b>Appendix C.</b> Examples of Completed Analyses .....	35
<b>The Quality Management System Approach</b> .....	51

Currently in preview, click buy full version

## Foreword

The primary goal of manufacturers of *in vitro* diagnostic (IVD) measurement procedures is to accurately report measurand values for patient samples. These measurement procedures may not produce accurate results when processed samples such as external quality assessment (EQA) samples, proficiency testing (PT) samples, or QC samples are measured. Because such processed sample matrixes typically undergo supplementation with additional components, and therefore are altered in some manner, measured results may not reflect the accuracy that would be observed for patient samples. Processed samples that have the same measurement response as that of patient samples that contain the same amount of an analyte are called commutable, while those that do not are called noncommutable. In this document, a matrix effect is defined broadly as differing test result biases in processed samples vs unprocessed patient samples due to unknown causes. The matrix effects that cause biases compared with unprocessed patient samples could be correlated to differences in conditions as encompassing as the entire measurement system or as specific as a reagent lot within a single IVD device.

Biases due to matrix effects in processed samples have the potential to affect the quality of patient care because of an inaccurate estimation of the accuracy of a measurement procedure. Depending on the intended use of the processed sample, the impact can range from negligible to serious. For example, a specific bias in a measuring interval verification sample set may have a different impact on the quality of patient care than the same bias in a QC sample. A measuring interval sample set matrix-related bias can directly affect the measuring interval allowed in patient sample results, whereas a QC matrix-related bias may affect the interpretation of QC results following a reagent lot change.

The objective of EP14 is to provide methods for identifying noncommutability so that improvements in measurement procedure analytical specificity and matrix compatibility may be considered. For example, a beneficial outcome of an evaluation may be a change in the processed sample's matrix or its additives, with an improvement in sample commutability. The techniques described in this guideline are also valid for testing the commutability of other samples, such as patient samples that have been processed (eg, added preservatives or spiking material, diluted, depleted, or frozen). Such samples, designed to substitute for unprocessed patient samples, are referred to as surrogate samples in CLSI document EP39.<sup>1</sup> EP14 will be helpful in explaining differences in test material results between measurement procedures, especially when such material serves as a basis for determining measurement procedure performance.

### Overview of Changes

This guideline was revised in 2022 under the Limited Revision Process and replaces the third edition of the guideline, which was published in 2014. Several changes were made in this edition, including:

- Adding content focused on processed samples such as EQA samples, PT samples, QC samples, and altered patient samples
- Adding an option for using preset commutability assessment criteria instead of prediction intervals
- Clarifying that EP14 is not designed for use with IVD manufacturer-specific calibrators and that the assessment methods described in this guideline should not be used in product regulatory submissions from such manufacturers

#### KEY WORDS

bias

commutability

Deming regression

interference

matrix

matrix effect

This page is intentionally left blank.

# Chapter ①

## Introduction

# Evaluation of Commutability of Processed Samples

## 1 Introduction

This section includes:

- Document scope and applicable exclusions
- Background information pertinent to the document content
- Standard Precautions information, as applicable
- Terms and definitions used in the document
- “Note on Terminology” that highlights particular use and/or variation in use of terms and/or definitions, where applicable
- Abbreviations and acronyms used in the document

### 1.1 Scope

This guideline provides protocols for evaluating commutability of processed samples when tested with quantitative measurement procedures. Such processed samples may be those created for proficiency testing/external quality assessment (PT/EQA), measuring interval verification sample sets, or QC samples. Processed samples can also be human specimens that are modified in a way that may change their measurement characteristics. In such cases, only a few processed samples (eg, three to five) should be evaluated to represent the behavior of the modification process being assessed. This restriction is suggested because the underlying analyses in this guideline do not account for the simultaneous assessment of many processed samples. CLSI document EP39<sup>1</sup> mentions the use of surrogate samples for measurement procedure performance testing, sometimes using many such samples, eg, for measurement procedure comparison studies. EP14 does not provide methods to screen a large number of samples in such cases.

This guideline is intended for developers of commercial diagnostic tests as well as laboratory-developed tests. This guideline is also useful for manufacturers of measuring interval sample sets and QC samples, and for PT or EQA providers to determine whether their processed samples are commutable with patient samples when examined with specified measurement procedures. Other options for assessing the commutability of processed samples across multiple measurement procedures are covered in CLSI document EP30<sup>2</sup> and in the literature.<sup>3</sup>

This guideline may also be useful to all clinical laboratory professionals wishing to investigate the commutability of a processed sample analyzed with a specific *in vitro* diagnostic (IVD) device in their laboratory.

EP14 is intended to assist in the education of clinical laboratorians and diagnostic manufacturers about the commutability of processed materials and how a sample’s matrix can affect some measurand values and their interpretation (referred to as matrix effects). For example, professionals may not be warned of a matrix effect caused by the interaction of processed PT/EQA material and the measurement procedure, and therefore the data may suggest to them that erroneous patient results are being generated, when in fact the results may be acceptable. Examples of a matrix effect due to the interaction of a processed QC and certain reagent lot(s) exist in the literature.<sup>4</sup> Therefore, these types of effects should not be a surprise to experienced laboratory staff and should not lead to erroneous conclusions about the suitability of results for patient samples. This guideline should assist all interested parties in not only evaluating the presence or absence of a matrix effect, but also by increasing awareness that the intended use of a processed matrix potentially affects patient care.