



CLINICAL AND  
LABORATORY  
STANDARDS  
INSTITUTE

2nd Edition

---

# CLSI EP32™

## Implementation of Metrological Traceability in Laboratory Medicine

CLSI EP32 provides guidance on establishing, validating, and documenting metrological traceability for end-user calibrators and results for human samples measured using *in vitro* diagnostic medical devices (IVD MDs) in medical laboratories based on the metrological traceability requirements for IVD MDs in ISO 17511.

---

A guideline for global application developed through the Clinical and Laboratory Standards Institute consensus process.

# Clinical and Laboratory Standards Institute

*Setting the standard for quality in medical laboratory testing around the world.*

The Clinical and Laboratory Standards Institute (CLSI) is a not-for-profit membership organization that brings together the varied perspectives and expertise of the worldwide laboratory community for the advancement of a common cause: to foster excellence in laboratory medicine by developing and implementing medical laboratory standards and guidelines that help laboratories fulfill their responsibilities with efficiency, effectiveness, and global applicability.

## Consensus Process

Consensus—the substantial agreement by materially affected, competent, and interested parties—is core to the development of all CLSI documents. It does not always connote unanimous agreement but does mean that the participants in the development of a consensus document have considered and resolved all relevant objections and accept the resulting agreement.

## Commenting on Documents

CLSI documents undergo periodic evaluation and modification to keep pace with advances in technologies, procedures, methods, and protocols affecting the laboratory or health care.

CLSI's consensus process depends on experts who volunteer to serve as contributing authors and/or as participants in the reviewing and commenting process. At the end of each comment period, the committee that developed the document is obligated to review all comments, respond in writing to all substantive comments, and revise the draft document as appropriate.

Comments on published CLSI documents are equally essential and may be submitted by anyone, at any time, on any document. All comments are managed according to the consensus process by a committee of experts.

## Appeal Process

When it is believed that an objection has not been adequately considered and responded to, the process for appeal, documented in the *CLSI Standards Development Policies and Processes*, is followed.

All comments and responses submitted on draft and published documents are retained on file at CLSI and are available upon request.

## Get Involved—Volunteer!

Do you use CLSI documents in your workplace? Do you see room for improvement? Would you like to get involved in the revision process? Or maybe you see a need to develop a new document for an emerging technology? CLSI wants to hear from you. We are always looking for volunteers. By donating your time and talents to improve the standards that affect your own work, you will play an active role in improving public health across the globe.

For additional information on committee participation or to submit comments, contact CLSI.

Clinical and Laboratory Standards Institute

P: +1.610.688.0100

F: +1.610.688.0700

[www.clsi.org](http://www.clsi.org)

[standard@clsi.org](mailto:standard@clsi.org)

# Implementation of Metrological Traceability in Laboratory Medicine

Jesper V. Johansen, PhD, MS  
Hubert W. Vesper, PhD  
Tony Badrick  
Ashley Beasley Green, PhD  
Jeffrey R. Budd, PhD  
Uliana Danilenko, PhD  
Denis Grote-Koska

Kevin Horan, BsC  
Anja Kessler  
W. Gregory Miller, PhD  
Nkemakonam C. Okoye, PhD, DABCC  
Mauro Panteghini, MD  
Elvar Theodorsson, MD, PhD

## Abstract

Clinical and Laboratory Standards Institute EP32—*Implementation of Metrological Traceability in Laboratory Medicine* provides guidance on establishing, validating, and documenting metrological traceability for end-user calibrators and results for human samples measured using *in vitro* diagnostic medical devices (IVD MDs) in medical laboratories based on the metrological traceability requirements for IVD MDs in ISO 17511.<sup>1</sup> Though CLSI EP32 is intended for use primarily by manufacturers of IVD MDs, the concepts and approaches recommended may be extended to apply to measurements performed in the medical laboratory either with commercially available or laboratory-developed tests.<sup>1</sup>

Clinical and Laboratory Standards Institute (CLSI). *Implementation of Metrological Traceability in Laboratory Medicine*. 2nd ed. CLSI guideline EP32. (ISBN 978-1-68440-268-7 [Print]; 978-1-68440-269-4 [Electronic]). Clinical and Laboratory Standards Institute, USA, 2025.

The Clinical and Laboratory Standards Institute consensus process, which is the mechanism for moving a document through two or more levels of review by the health care community, is an ongoing process. Users should expect revised editions of any given document because rapid changes in technology may affect the procedures, methods, and protocols in a standard or guideline, and users should replace outdated editions with the current editions of CLSI documents. Current editions are listed in the CLSI catalog and posted on our website at [www.clsi.org](http://www.clsi.org).

**If you or your organization is not a member and would like to become one, or to request a copy of the catalog, contact us at:**

**P:** +1.610.688.0100 **F:** +1.610.688.0700 **E:** [customerservice@clsi.org](mailto:customerservice@clsi.org) **W:** [www.clsi.org](http://www.clsi.org)

Copyright ©2025 Clinical and Laboratory Standards Institute. Except as stated below, any reproduction of content from a CLSI copyrighted standard, guideline, or other product or material requires express written consent from CLSI. All rights reserved. Interested parties may send permission requests to [permissions@clsi.org](mailto:permissions@clsi.org).

CLSI hereby grants permission to each individual member or purchaser to make a single reproduction of this publication for use in its laboratory procedures manual at a single site. To request permission to use this publication in any other manner, e-mail [permissions@clsi.org](mailto:permissions@clsi.org).

To read CLSI's full Copyright Policy, please visit our website at <https://clsi.org/terms-of-use/>.

## Suggested Citation

CLSI. *Implementation of Metrological Traceability in Laboratory Medicine*. 2nd ed. CLSI guideline EP32. Clinical and Laboratory Standards Institute; 2025.

### Previous Edition:

February 2006

CLSI EP32-Ed2

ISBN 978-1-3440-268-7 (Print)

ISBN 978-1-68440-269-4 (Electronic)

ISSN 1558-6502 (Print)

ISSN 2162-2914 (Electronic)

Volume 45, Number 6

## Committee Membership

### Consensus Council

The Consensus Council sets priorities for CLSI standards development and votes on Final Draft documents to confirm that process requirements have been met. Consensus Council members are listed on the CLSI website: <https://clsi.org/standards-development/consensus-council/>

### Document Development Committee on Metrological Traceability

<b>Jesper V. Johansen, PhD, MS</b> <b>Chairholder</b> <b>Radiometer Medical ApS</b> <b>Denmark</b>	Thomas Ciesiolka, PhD Roche Diagnostics GmbH Germany	Majid Moridani Laboratory Corporation of America USA
<b>Hubert W. Vesper, PhD</b> <b>Vice-Chairholder</b> <b>Centers for Disease Control and Prevention</b> <b>USA</b>	Ashley Beasley Green, PhD National Institute of Standards and Technology USA	Mauro Pantellini, MD Department of Laboratory Medicine, Ludwik Rydygier Collegium Medicum in Bydgoszcz, Nicolaus Copernicus University Poland
Tony Badrick RCPA Quality Assurance Programs Pty Limited Australia	Anja Kessler Reference Institute for Bioanalytics (RfB) Germany	Elvar Theodorsson, MD, PhD Linköping University Sweden

### Expert Panel on Evaluation Protocols

Expert panel volunteers support the development of CLSI documents by providing technical expertise in specialty areas. Expert panel members are listed by area of expertise on the CLSI website: <https://clsi.org/standards-development/expert-panels/>

## Acknowledgment

---

CLSI, the Consensus Council, and the Document Development Committee on Metrological Traceability gratefully acknowledge the following volunteers for their important contributions to the revision of CLSI EP32:

Jeffrey R. Budd, PhD  
USA

Denis Grote-Koska  
Medizinische Hochschule Hannover/  
RfB  
Germany

W. Gregory Miller, PhD  
Virginia Commonwealth University  
Health System  
USA

Uliana Danilenko, PhD  
Centers for Disease Control and  
Prevention  
USA

Kevin Horan, BSc  
Siemens Healthcare  
USA

Nkemakonam C. Okoye, PhD, DABCC  
Northwell Health Laboratories  
USA

# Contents

Abstract	i
Committee Membership	iii
Foreword	vii
<b>Chapter 1: Introduction</b>	<b>1</b>
1.1 Scope	1
1.2 Background	2
1.3 Standard Precautions	4
1.4 Terminology	4
<b>Chapter 2: Path of Workflow</b>	<b>11</b>
<b>Chapter 3: Metrological Traceability in Laboratory Medicine</b>	<b>13</b>
<b>Chapter 4: Calibration Hierarchies</b>	<b>15</b>
4.1 A Model Calibration Hierarchy for an End-User Measurement Procedure	16
4.2 A Realistic Calibration Hierarchy for End-User Measurement Procedures in Laboratory Medicine	18
4.3 References in Metrological Traceability	21
4.4 Certified Reference Materials	22
4.5 Reference Measurement Procedures	26
4.6 Measurement Uncertainty	27
<b>Chapter 5: Commutability</b>	<b>29</b>
5.1 Explanation of Which Calibrators Need to be Commutable With Human Samples	31
5.2 Handling Noncommutability of Reference Materials	31
<b>Chapter 6: Requirements for Establishing Metrological Traceability</b>	<b>37</b>
6.1 Introduction	38
6.2 Definition of the Measurement Standard	38
6.3 Specifications for Maximum Allowable Expanded Measurement Uncertainty	39
6.4 Defining the Calibration Hierarchy	39
6.5 Selection and Requirements for Reference Materials and Calibrators	50
6.6 Selection and Requirements for Measurement Procedures	50
6.7 Estimating Measurement Uncertainty of Assigned Values for <i>In Vitro</i> Diagnostic Medical Devices	51
6.8 Validation of Metrological Traceability of Values Assigned to an End-User <i>In Vitro</i> Diagnostic Medical Device Calibrator and to Human Samples	52
6.9 Additional Calibration Hierarchy Documentation Responsibilities	53

## Contents (Continued)

Chapter 7: Conclusion .....	55
Chapter 8: Supplemental Information .....	57
References .....	58
Additional Resources .....	64
Appendix A. Overview of Worked Examples .....	65
Appendix B. Worked Example of Glucose in Plasma .....	68
Appendix C. Worked Example of pH in Whole Blood .....	77
Appendix D. Worked Example of Alanine Aminotransferase Catalytic Concentration in Serum .....	105
Appendix E. Worked Example of Free Thyroxine in Serum .....	110
Appendix F. Worked Example of Hemoglobin A <sub>1c</sub> in Whole Blood .....	116
Appendix G. Worked Example of Human Chorionic Gonadotropin and Immunoglobulin G in Serum .....	124
Appendix H. Worked Example of D-dimer in Serum .....	132
The Quality Management System Approach .....	146

## Foreword

---

CLSI EP32 provides guidance for *in vitro* diagnostic medical device (IVD MD) developers to implement metrological traceability in laboratory medicine and describes the benefits associated with its implementation.

CLSI EP32 aids in interpreting and implementing metrological traceability according to ISO 17511.<sup>1</sup> It also references the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC) recommendations on assessing commutability of calibrators, when applicable, because their commutability is crucial for implementation of calibration hierarchies. CLSI EP32 covers the principles of all six model calibration hierarchies described in ISO 17511<sup>1</sup> and their implementation in situations representing substantial variations of the six model calibration hierarchies.

CLSI EP32 outlines what is required by IVD MD developers to establish metrological traceability. It provides guidance on explaining the results of studies to end users and describes what end-user laboratories must do to verify results based on traceability concepts. These concepts are explained in general terms in the main text and in detailed examples in the appendixes.

Establishing, validating, and documenting metrological traceability is part of a range of activities for improving and maintaining the usefulness of measurements performed in laboratory medicine and public health. Implementing metrological traceability is facilitated by different frameworks, such as those developed by formal national and international standardization programs or stakeholder organizations.

Examples of such frameworks are those from the Centers for Disease Control and Prevention's standardization programs for lipids, hormones, and vitamin D<sup>2</sup>; the NGSP for hemoglobin A<sub>1c</sub> (HbA<sub>1c</sub>),<sup>3</sup> and the IFCC (HbA<sub>1c</sub><sup>4</sup> and thyroxine).<sup>5</sup>

Reference measurement laboratories operating reference measurement procedures to assign values to human samples (HS) or calibration materials and trueness controls must provide information about the highest order of reference to which their measurement result is traceable and how this traceability was established, according to ISO 15195.<sup>6</sup>

### Overview of Changes

CLSI EP32-Ed2 replaces the previous edition of the approved report, CLSI EP32-R, published in 2006. CLSI EP32 has been completely rewritten and focuses on reflecting the changes implemented in the current edition of ISO 17511,<sup>1</sup> specifically:

- Metrological traceability of measurement results for HS, not just values assigned to product calibrators
- Details on which calibrators in a calibration hierarchy must be commutable and how to handle noncommutability in these materials
- Discussion of validation of metrological traceability

Furthermore, seven worked examples have been developed to illustrate key aspects and considerations for different

calibration hierarchies. An introduction to the worked examples is found in Appendix A. The worked examples are included as Appendixes B through H:

- Appendix B: Worked Example of Glucose in Plasma
- Appendix C: Worked Example of pH in Whole Blood
- Appendix D: Worked Example of Alanine Aminotransferase Catalytic Concentration in Serum
- Appendix E: Worked Example of Free Thyroxine in Serum
- Appendix F: Worked Example of Hemoglobin A<sub>1c</sub> in Whole Blood
- Appendix G: Worked Example of Human Chorionic Gonadotropin and Immunoglobulin G in Serum
- Appendix H: Worked Example of D-dimer in Serum

**NOTE:** The content of CLSI EP32 is supported by the CLSI consensus process and does not necessarily reflect the views of any single individual or organization.

**KEY WORDS**

calibration hierarchy

calibrator

certified reference material

commutability

measurement uncertainty

metrological traceability

reference measurement procedure

traceability chain

value assignment

# Chapter ①

## Introduction

# Implementation of Metrological Traceability in Laboratory Medicine

---

## 1 Introduction

### 1.1 Scope

CLSI EP32 provides information for implementing metrological traceability according to ISO 17511.<sup>1</sup> It describes the necessary components and their use for *in vitro* diagnostic medical device (IVD MD) manufacturers to correctly establish, implement, and maintain metrological traceability. CLSI EP32 explains different calibration hierarchies for end-user medical laboratory measuring systems that are metrologically traceable to the highest available measurement procedures (MPs) and calibration materials.

CLSI EP32 describes related procedures (eg, commutability assessment, estimation of measurement uncertainty [MU], and analytical performance specifications [APS]) used in implementing metrological traceability. However, CLSI EP32 refers to other publications for detailed design and execution of these procedures.

### 1.2 Background

A primary goal of laboratory medicine is to provide test results to support clinical decision-making and foster optimal health care. The test results should be interpretable to end users regardless of the laboratory or IVD MD. The test results must be equivalent for the same measurand from various MPs and laboratories.

Metrological traceability is necessary for the establishment and maintenance of the equivalence of measurement results. Metrological traceability is a property of a measurement result that can be related to a higher-order reference material (RM) or reference measurement procedure (RMP) through a documented unbroken chain of calibrations, each contributing to the MU.

A reference selected for metrological traceability should have a characteristic(s) among the following:

- Link of measured quantities to the definition of the International System of Units (SI)
- A certified value of an RM
- The value assigned using an RMP
- The value assigned to an international conventional calibrator (ICC)
- The values assigned through an international harmonization protocol

If none of the above characteristics is available, an IVD MD developer is responsible for defining its internal reference.

Measurements of distances, masses, and time are used routinely in many aspects of daily life and industry. Therefore, such measurements are often intuitively understood. In contrast, measurements in laboratory medicine are usually concerned with the concentration of molecules frequently present in very low concentrations among a multitude of other molecules that comprise the human sample (HS) matrix. Furthermore, the molecules are rarely measured directly. Instead, measurement signals related to the analyte concentrations in HS are created through selective chemical reactions to target the molecules of interest. These reactions can be detected using physical methods, commonly the absorbance or emission of light at specific wavelengths, voltages, or an electric current as illustrated in Figure 1.